RapidPlan

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
RAPIDPLAN™—KNOWLEDGE-BASED PLANNING: PRODUCT RATIONALE

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
PQM Score Verteilung


* Variation in external beam treatment plan quality: An inter-institutional study of planners and planning systems. Nelms, et al; Practical Radiation Oncology, Vol2, Issue 4, Oct.-Dec 2012

KNOWLEDGE-BASED PLANNING*: IMPLEMENTATION

* Data using knowledge based planning software similar to RapidPlan. Moore, K.L. *et al.* Experience-Based Quality Control of Clinical Intensity-Modulated Radiotherapy Planning. IJROBP 81(2): 545-551 (2011)

KNOWLEDGE-BASED PLANNING*: IMPLEMENTATION

* Data using knowledge based planning software similar to RapidPlan. Moore, K.L. *et al.* Experience-Based Quality Control of Clinical Intensity-Modulated Radiotherapy Planning. IJROBP 81(2): 545-551 (2011)

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eoc

Luca Cozzi
Antonella Fogliata

IOSI, Bellinzona,
Switzerland

RapidPlan model for SIB Head and Neck


- Very inhomogenous dataset
- SIB head and neck with different dose prescriptions and different number of PTVs (2 or 3)
 - 16 pts with 3 PTVs: H) 65.4-71.9Gy; M) 64-54Gy; L) 50-59Gy
 - 24 pts with 2 PTVs: H) 63.6-69.6Gy; L) 52-62Gy
- 40 patients included in the model
- Almost all plans to configure the model were IMRT
- Structures in the model: 2 or 3 PTVs, spinal cord, parotid left, parotid right
- All organs at risk had generated priorities: line and mean dose for parotids, max dose for spinal cord

RAPIDPLAN ACTIVITIES OF IOSI BELLINZONA

Model based vs. original plan

▲ Model based RapidArc (69.96Gy)

■ Original IMRT plan (66.6Gy)

Mean right parotid: 25.2 Gy
Mean left parotid: 17.4 Gy
Max sp.cord: 34.1 Gy

Mean right parotid: 41.7 Gy
Mean left parotid: 23.8 Gy
Max sp.cord: 34.1 Gy

MODEL BASED PLAN BETTER THAN ORIGINAL

RAPIDPLAN™ CLINICAL COLLABORATORS FOCUS GROUP 2014: RAPIDPLAN: EVALUATING CONSISTENCY BETWEEN CAMPUSES USING AN ESOPHAGUS MODEL

Sean L Barry
Amanda Boczkowski Rongtao Ma
Penpeng Zhang
Margie Hunt

Memorial Sloan-Kettering
Cancer Center (New York)

Motivation


- Want to use knowledge based planning (KBP) to identify and reduce treatment plan variation within an institution


RAPIDPLAN: EVALUATING CONSISTENCY BETWEEN CAMPUSES USING AN ESOPHAGUS MODEL

Building the model

- Data Source:
 - > Retrospective selection of *main campus* IMRT esophagus plans created b/w 2009-2014.
 - > Started w/ 64 candidate plans
 - Extracted dosimetric details and reviewed/excluded:

Non-standard anatomy

Non-standard clinical trade-offs

- > Ultimately, 58 patients were available for the RapidPlan Model Configuration

RAPIDPLAN: EVALUATING CONSISTENCY BETWEEN CAMPUSES USING AN ESOPHAGUS MODEL

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Memorial Sloan Kettering
Cancer Center.

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Oncology Institute of Southern Switzerland (Bellinzona, Switzerland) and Cheng-Ching General Hospital (Taichung, Taiwan).

Aim of the study

Evaluation of a knowledge-based dose-volume constraints prediction engine recently implemented in a commercial treatment planning system (TPS)

Site Selection


Advanced stage hepatocellular carcinoma (HCC) was chosen due to high variability in tumor size and position in the liver.

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RAPIDPLAN™ — ENTHALTENE MODELLE

Cancer Care Manitoba

- CCMB Head & Neck Model
- CCMB Prostate Model
- HFHS Lung SBRT Model
- WUSTL Head & Neck Model
- WUSTL Prostate Lymph Node Model
- WUSTL Prostate Model

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
Washington University

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RAPIDPLAN™—KNOWLEDGE-BASED PLANNING BENEFITS

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RAPIDPLAN™—KNOWLEDGE-BASED PLANNING

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